

High-performance liquid chromatography

Detectors

fall 2015

1

Requirements for detectors (1)

- Equally sensitive towards all components.
OR
Sensitive towards components of interest.
- Gives as much information about the analyte as possible.
- Independent of eluent composition (gradient) and temperature.
- Detects very low concentrations (analysis of traces).
- Does not cause widening of the peaks (small flow cell).

fall 2015

2

Requirements for detectors (2)

- Fast, in order to “catch” very narrow and fast moving peaks.
- Stable and reproducible signal.
- Better if wide linear range (at least 3 orders of magnitude).
- Non-destructive.
- Easy to use, robust.
- Cheap.

fall 2015

3

Types of detectors (1)

- Concentration and mass selective:
 - **Concentration selective** – signal is dependent of the analyte concentration in eluate.
 $S \propto c \text{ (g ml}^{-1}\text{)}$
 - **Mass selective** – signal is dependent on the mass flow, amount of sample molecules in time.
 $S \propto n/\Delta t \text{ (g s}^{-1}\text{)}$
- Eg when switching of the pump at the highest point in chromatographic peak – for concentration selective detector signal remains, for mass selective detector signal will fall to the base noise level.

fall 2015

4

Types of detectors (2)


- Selectivity of the measured properties:
 - Non-selective detectors react to the **bulk property of the passing solution**
 - Refractive index
 - Conductivity
 - Selective detectors measure a response due to specific property of the **analyte molecule**
 - UV absorbance
 - Fluorescence
 - Mass to charge ratio

fall 2015

5

Sensitivity


For the estimation of sensitivity it is not enough to evaluate the signal of the analyte. Sensitivity is estimated through the signal to noise ration.



fall 2015

Flow cell (1)


- Cells with 1-15 μl volume are used:
 - Commonly used HPLC flow cells are 6-10 μl
 - UHPLC micro-flow cells are maximum 2.7 μl



fall 2015

7

Flow cell (2)



fall 2015

8


UV-Vis detector

- Most widely used HPLC detector.
 - Selective detector (measures properties of the analyte).
 - Concentration selective.
 - Sensitive.
 - Wide linear range (4 orders of magnitude).
 - Not sensitive to temperature changes.
 - Suitable for gradient elution.

fall 2015

9

Electromagnetic spectrum




UV Cut-off - for common organic solvents used in reverse phase HPLC:
Acetonitrile 210nm / Methanol 210nm / THF 220nm

fall 2

10

UV-Vis: principle

- Sometimes double-beam systems with reference light passing through the reference cell and filter are used.



$$A = \log \frac{I_0}{I}$$

- Beer's law: $A = \epsilon l c$ ehk $A = abc$, where A – absorbance, ϵ (a) – molar absorptivity, l (b) – optical path length, c – analyte concentration.

fall 2015

11

UV-Vis: types and radiation sources

Fixed wavelength and filters (1)


- In case of low pressure Hg lamps wavelength 254 nm is used.
- Cd-lamp – 229 nm.
- Zn-lamp – 214 nm.
 - Other wavelengths are filtrated out. (Why?)
 - Significantly more sensitive than detectors with changeable wavelength.

fall 2015

12

UV-Vis: types and radiation sources Fixed wavelength and filters (2)

- Deuterium-lamp – constant UV-specter 190 nm to ca 340 nm.



fall 2015


13

UV-Vis: types and radiation sources Fixed wavelength and filters (3)

- Tungsten lamp – 340-850 nm.

Lamp bulb is made of quartz and can be used at temperatures as high as 900 C and pressure at 20 bars.

To lengthen the lifetime of lamp, the bulb is filled with a halogen such as Br or I.



fall 2015

14

UV-Vis: applications


- Used for detection of very many compounds:
 - Double bond and free electron pair atom closeby $X=Y-Z$
 - Br-, I- and S-containing compounds.
 - Carbonyl (C=O) or nitro (NO_2).
 - Conjugated double bond $X=X-X=X$.
 - Aromatic ring.
 - Anions: Br^- , I^- , NO_3^- , NO_2^-
- **Not suitable** for detection of saturated hydrocarbons and their amino and nitrile derivatives..

fall 2015

15

UV-Vis: types and radiation sources

- With monochromator (D_2 - and/or W-lamp)
- Diode array - PDA-photodiode array, DAD-diode array detector



fall 2015

16


UV-Vis: working modes

- Detection of chromatogram at fixed wavelength – all listed types.
- Changing the wavelength during detection of the chromatogram – most detectors with D_2 - and W- lamps.
- Scanning (D_2 and W) detectors allow to stop the eluent flow and then register the spectrum.
- Diode array allows to detect **full spectrum** at each point in the chromatogram (with eluent flow).
 - In addition to sample wavelength, reference wavelength can be determined.

fall 2015

17

UV-Vis: choosing the wavelength




fall 2015

18

UV-Vis: Reference wavelength

- For diode array detector a reference wavelength (range) can be determined that corrects the drift for example from gradient.




fall 2015

19

FLD – fluorescence detector (1)

- FLD – excitation radiance is perpendicular to the detection:



fall 2015

Hitachi High-Tech

FLD – fluorescence detector (2)


- Properties
 - Selective detector.
 - LOD is **1000 times** better than in UV-Vis.
 - For higher sensitivity, flow cell volume is *ca* 20 μ l.
 - Sensitive to **fluorescence depleting** components in the mobile phase (eg dissolved oxygen).
 - Linear range depends of the system, but is narrow.

fall 2015

21

RI – refractive index detector (1)

- RI – *refractive index*



fall 2015

Hitachi High-Tech

RI – refractive index detector (2)

- Properties
 - Non-destructive detector (refractive index of the whole eluate).
 - LOD is 1000 times worse than for UV-Vis detector.
 - Very sensitive towards temperature changes – needs a thermostate.
 - It is necessary to have a reference cell with an eluent.
 - Not suitable for gradient elution.
 - Gases dissolved in eluent will influence the detection.

fall 2015

23

ELSD – evaporative light scattering detector


- ELSD – eluent with the dissolved analytes are nebulized and evaporated into tiny sample particles. These scatter a light beam (laser) and the extent of scattering is measured.
 - Non selective detector (sensitive towards **all non volatile** compounds).
 - More sensitive than RI.
 - Can be used for gradient elution if eluent contains volatile components.
 - Linear range is basically inexistent.

fall 2015

24

IR – infrared absorbance detector

- Based on the property of organic compounds to absorb infrared light.
 - Selective detector.
 - Mobile phase must be chosen so that it does not absorb at the same wavelength as the analyte.
 - Volatile buffers.
 - Good for polymer analysis.
 - Ability to differentiate secondary structure.



fall 2015

Electrochemical detectors


- Used for different types:
 - Amperometric
 - Polarographic
 - Coulometric
 - Conductometric

fall 2015

26

Amperometric (1)

- Amperometric detector measures electrical current resulting from oxidation or reduction reactions.



fall 2015



Amperometric (2)

- Properties
 - Selective detector (as low as fmole/L).
 - LOD is very low.
 - Simple and cheap.
 - Wide linear dynamic range (more than 6 orders of magnitude)
 - Mobile phase must conduct electricity, but does not have to be aqueous (non-polar eluents are not suitable).
 - Sometimes difficult to find the best conditions.

fall 2015

28

Conductivity detector (1)

- Classical detector for ion chromatography.
- Conductivity of the mobile phase is measured.



fall 2015

29

Conductivity detector (2)

- Properties
 - Non-selective mass sensitive detector.
 - Signal is dependent on of the ionic component in the mobile phase.
 - **Sensitivity is decreased** if eluent is conducting electricity (suppression columns are used).
 - Sensitive to temperature.
 - Linear range is not very wide.

fall 2015

30

MSD – mass spectrometric detector

- MS as a detector is so powerful and special that the whole method is abbreviated as LC-MS (LC/MS).
- MS has developed so much in the last decade that options presented in older books must be taken with care.
- For MS detection, analyte must be in ionic form, therefore MS is splitted into to parts:
 - Ion source
 - Mass spectrometer

fall 2015

31

MSD properties

- **Selectivity** – Coeluting peaks can be separated/isolated according to mass, even when they are not chromatographically separated.
- **Detection of peaks** – specific to chemical compounds.
- **Information on molecular mass** – Identification of known and unknown compounds.
- **Information on the structure of the compound** – Fragmentation provides additional information.
- **Faster method development** – Fast and easy analyte detection without retention time validation.
- **Differentiating sample matrix from analyte** – may reduce the time spent on sample preparation.
- **Quantification** – Both quantitative and qualitative analysis.

fall 2015

32

MSD ion sources

- Roles of ionization source:
 - Ionization of the analyte.
 - Separation of analyte from the eluent.
 - Make sure that only analyte enters the mass spectrometer.
- Atmospheric ionization sources (API):
 - ESI (*electrospray*)
 - APCI (*AP chemical ionization*)
 - APPI (*AP photoionization*)
 - AP-MALDI (*AP matrix-assisted laser desorption/ionization*)

fall 2015

33

MSD: ESI

- Principle of ESI is electrospray – eluent sprayed into high voltage creates an aerosol. Small drops in the aerosol will acquire a charge.
- Eluent flow introduced to ESI is sprayed into fine droplets with gaseous N₂ and electric field. Analyte molecules acquire a charge (eg addition of H⁺) and due to electric field move to MS.
- ESI
 - **Very soft** ionization, even biomolecules don't decompose.
 - Most commonly used ionization source.
 - Concentration sensitive.

fall 2015

34

MSD: APCI

- In case of APCI eluent flow is sprayed to droplets with the help from gaseous N₂. Eluent is evaporated through heating.
- Corona discharge is used to ionize N₂, which will ionize eluent molecules and which in turn will ionize analyte.
- APCI
 - **Very soft** ionization, which allows to analyze **less polar compounds** than ESI.
 - Most common after ESI.
 - Can be said that mass sensitive.

fall 2015

35

MSD: mass spectrometers

- Quadrupole MS.
 - Single quadrupole is quite sensitive with good LOD but it does not allow fragmentation and is therefore less selective than other MS types.
 - Triple quadrupole (3Q) is more sensitive than single quadrupole. MS² allows to get **structural information** and therefore raises the selectivity.
- Ion trap MS
 - Very sensitive MS detector; allows at least MS⁴, which adds significant amount of **structural information**.
- Time of flight (TOF) MS
 - Very exact mass.
- Orbitrap MS
 - Most sensitive MS detector, very exact for small and large molecules. Fast detection, wide linear range and allows MSⁿ which gives a lot of **structural information**.

fall 2015

36


More detectors

- Photoconductivity detector
- Radioactivity detector
- Light scattering detector
- Dielectric constant detector
- Corona Charged Aerosol Detector: http://www.esainc.com/products/type/hplc_systems/detectors/coronacad

fall 2015

37

CCAD (1)



fall 2015

38

CCAD (2)

- Excellent sensitivity (in nanograms)
- Wide dynamic range
- Very reproducible signal
- Can detect analytes with very different properties (also without chromophores)
- Robust and easy to use

fall 2015

39